# **Conversion to BIDS Format**

Converting data to the BIDS (Brain Imaging Data Structure) format involves several steps to properly organize and structure your data. Below is a clear step-by-step guide:

# 1. Prepare the Initial Folder Structure

- Create a main folder called DATA.
- Inside DATA, create two folders: bids and raw.
- Place your DICOM images inside a folder named (for example) dcm within the raw folder.

### 2. Navigate to the Directory

• Go to the directory where your raw folder is located.

cd raw

#### 3. Check the Contents of the raw Folder

List the contents of the raw folder.
 You should see subfolders named like 01, 02, 03, etc., each corresponding to a different subject.

dir
(or ls on Unix-based systems)

#### 4. Install DicomSorter

Install the dicomsorter package using pip:

pip install dicomsorter

# 5. Sort Directories Using DicomSorter

- Use dicomsorter to organize your DICOM files.
- Replace . /01 / with the path to your subject's DICOM folder, and . /sub-01 / with the desired output folder.

dicomsorter ./01/ ./sub-01/

#### 6. Install BidsCoin

• Install the bidscoin package using pip:

pip install bidscoin

#### 7. Go Back to the DATA Folder

• Return to the main DATA folder.

cd ..

# 8. Map Raw Data to BIDS Structure

- Use bidsmapper to map your raw data into the BIDS folder.
- Adjust your configuration files according to your processing needs.

bidsmapper raw bids

# 9. (Optional) Edit the Mapping

• Open the BIDS editor to review and modify the mapping if necessary.

# bidseditor -b ./bids/code/bidsmap bids

#### 10. Finalize the Conversion to BIDS

• Use bidscoiner to complete the conversion of your data to the BIDS format.

bidscoiner raw bids

# **Important Notes:**

- Replace all folder names and paths with the ones you are using.
  Pay attention to any error messages or additional requirements mentioned in the specific tutorials you follow.